

Performance Comparison: Ibsen Photonics PEEBLE NIR (PBM-400) vs. similar Compact Spectrometer

In this document, the findings from comparative measurements carried out with a PEEBLE NIR spectrometer and a comparable compact spectrometer for reference are summarized.

Device Specifications

The table below shows a comparison of some important specifications of each device:

Parameter	PEEBLE NIR	Reference Compact Spectrometer
Approx. Dimensions / mm	40 x 23 x 13	80 x 60 x 10
Approx. Weight / g	< 5	< 90
Number of Pixels	256 (128 are used)	256
Slit Width / μm	25	25
Detector Type	uncooled InGaAs	uncooled InGaAs
Integration Time / μs	48 – 100000	1 – 100000
Numerical Aperture	0.22	0.22
Fiber Connection	SMA905	SMA905
Spectral Range / nm	950 – 1700	950 – 1700
Spectral Resolution / nm	16	7
Operating Temperature / $^{\circ}\text{C}$	0 – 50	5 – 50
Storage Temperature / $^{\circ}\text{C}$	-20 – 70	-20 – 70

Experimental

For an experimental comparison, a simple setup as shown in Figure 1 was used (taken from Zimmerleiter et al. <https://publikationen.bibliothek.kit.edu/1000128686>, 2021).




Figure 1: Schematic drawing of the experimental setup used for the conducted comparison measurements.

For connecting both spectrometers, a round to linear fiber was used and carefully aligned to illuminate as much of the entrance slit as possible.

The exposure time was set to a level that was only slightly below saturation of the detector (with no sample present). Then averaging was adjusted to get a measurement duration (including storing the data on the PC hard drive) of approximately 2 seconds. The results for both devices are:

PEEBLE NIR:	Exposure time = 370 μ s Averaging = 2000
Reference Compact Spectrometer:	Exposure time = 800 μ s Averaging = 30

It seems that the PEEBLE NIR spectrometer outperforms the reference compact spectrometer when it comes to data transfer speed. This is most likely due to the possibility of averaging spectra directly in hardware, which was not possible with the other compact spectrometer.

100% Lines

To test for the noise level, the light straight from the halogen lamp was measured (no cuvette). 101 spectra measured straight from the halogen light source were acquired for each spectrometer. The first measured spectrum (S_0) was used as background and 100 individual 100%-lines were calculated with the formula below:

$$A_i^{100\%} = \log_{10}(S_0) - \log_{10}(S_i)$$

The 100%-lines for the two devices (200 lines in total) are shown in Figure 2. The left graph is a zoomed version of the one on the right side.

To compare the noise performance of the two spectrometers the root mean square was calculated according to:

$$\text{RMS} = \frac{1}{100} \sum_{i=1}^{100} \left[\frac{1}{N_\lambda} \sum_{j=1}^{N_\lambda} [A_i^{100\%}(\lambda_j)]^2 \right]^{1/2}$$

where N_λ is the total number of spectral points acquired.

The RMS-values of the two spectrometers are:

$$\begin{aligned} \text{RMS}_{\text{PEEBLE}} &= 1.03 \times 10^{-4} \\ \text{RMS}_{\text{COMP}} &= 2.28 \times 10^{-4} \end{aligned}$$

From the RMS values it can be seen, that the PEEBLE outperforms the reference compact spectrometer by approximately a factor of 2 when it comes to noise performance. However, it offers only approximately half the spectral resolution.




Figure 2: 100% Lines of the PEEBLE NIR compared to the reference compact spectrometer with similar specifications.

To make the spectra more comparable, a smoothing filter (Svaitzky-Golay, 3 points) was applied to the compact spectrometer data. This gave spectra with half the spectral resolution, which are comparable to the PEEBLE. The resulted RMS of approximately 2.20×10^{-4} was only slightly better.

A significant increase in performance is gained when using only the data between 1000 nm and 1650 nm, resulting in an RMS of 2.71×10^{-5} for the PEEBLE NIR and an RMS of 5.39×10^{-5} (4.62×10^{-5} with Savitzky-Golay filter) for the reference compact spectrometer.

Absorption Spectrum of Ethanol

A real-world comparison was done by acquiring 100 transmission spectra of ethanol in a quartz-glass cuvette with a pathlength of 1 mm with each instrument. The empty cuvette was used as a background. The resulting spectra (200 in total) are shown Figure 3.




Figure 3: Absorption of ethanol measured in a quartz-glass cuvette with 1 mm path length.

The acquired spectra agree nicely, with a slight offset of approximately 5 nm. However, due to the lower spectral resolution of the PEEBLE NIR some spectral features are not as clearly visible in the absorption spectra as with the reference compact spectrometer.

Comparison conducted and report written by Robert Zimmerleiter on March 2nd, 2022.

Contact data:

DI Robert Zimmerleiter

RECENDT – Research Center for Non-Destructive Testing GmbH
Altenberger Straße 69, Science Park 2
4040 Linz
AUSTRIA
Tel.: +43 732 2468 4624
Fax: +43 732 2468 4606
E-Mail: robert.zimmerleiter@recendt.at
<http://www.recendt.at>

DI Dr. Markus Brandstetter

RECENDT – Research Center for Non-Destructive Testing GmbH
Altenberger Straße 69, Science Park 2
4040 Linz
AUSTRIA
Tel.: +43 732 2468 4620
Fax: +43 732 2468 4606
E-Mail: markus.brandstetter@recendt.at
<http://www.recendt.at>

DI Dr. Christian Rankl

RECENDT – Research Center for Non-Destructive Testing GmbH
Altenberger Straße 69, Science Park 2
4040 Linz
AUSTRIA
Tel.: +43 732 2468 4644
Fax: +43 732 2468 4606
E-Mail: christian.rankl@recendt.at
<http://www.recendt.at>